

## Oral Fluid Drug Panels

<b>Oral Fluid 6-Panel (S6P)</b>	<b>Screening Cutoff (Dilute)</b>	<b>Confirmation Cutoff (Neat)</b>
Amphetamine <i>(Adderall, Benzedrine, Dexedrine, Vyvanse, etc.)</i>	100 ng/ml	50 ng/ml
Methamphetamines <i>Methamphetamine MDMA, MDA, MDEA</i>	40 ng/ml	50 ng/ml
Cocaine <i>Cocaine, Benzoylcegonine</i>	5 ng/ml	8 ng/ml
Opiates <i>Codeine Morphine Hydrocodone (Vicodin, Lortab/Lorcet) Hydromorphone (Dilaudid) 6-MAM (Heroin)</i>	10 ng/ml	40 ng/ml 4 ng/ml (6-MAM)
Phencyclidine (PCP)	1 ng/ml	10 ng/ml
Marijuana (THC)	1 ng/ml	2 ng/ml
<b>Oral Fluid 10-Panel (S10P) (In Addition to 6-Panel)</b>	<b>(Dilute)</b>	<b>(Neat)</b>
Oxycodone <i>Oxycontin, Percodan/Percocet</i>	10 ng/ml	40 ng/ml
Methadone	5 ng/ml	15 ng/ml
Benzodiazepines <i>Diazepam/Nordiazepam (Valium) Alprazolam/ Alpha-hydroxyalprazolam (Xanax)</i>	1 ng/ml	5 ng/ml
Barbiturates <i>Secobarbital (Seconal) Butalbital (Fioricet, Fiorinal) Pentobarbital (Nembutal) Phenobarbital (Luminal) Amobarbital (Amytal, Tuinal)</i>	20 ng/ml	60 ng/ml
<b>Separate Panels (Offered separate from 6 and 10 panels)</b>	<b>Screening</b>	<b>Confirmation (Neat)</b>
Synthetic Marijuana (SSTHC) <i>JWH-018, JWH-018NPA, JWH-018NHP, JWH-073, JWH-073NBA, JWH-073NHB, JWH-250, AM 2201</i>	1.0 ng/ml (Neat)	1.0 ng/ml
Bath Salts (SBS) <i>Mephedrone, MDPV, Methcathinone, Methylone</i>	1.0 ng/ml (Neat)	1.0 ng/ml
Cotinine (Nicotine Metabolite) (SCOT)	10.0 ng/ml (Dilute)	10.0 ng/ml

**Neat Vs. Dilute:**

Screening cutoff levels listed as Dilute to account for device buffer.  
Confirmation listed as Neat in accordance with SAMHSA guidelines.

**Validity Testing:**

All oral fluid panels include IgG human validity testing.

**Dose Response / Interpretation:**

Results should be considered as a positive or negative only. Levels of drugs found in oral fluid testing cannot be reliably used to make assumptions about donor usage.